

Emergency lower leg

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM03450



Model of a human leg with deep tissue injury, caused by an explosive (blind grenade).

Main advantages:

- The vessel system can be connected with a Luer-lock connector.
- The whole model is made of strong silicone material and can be used for tamponade treatment of the wounds.

Emergency leg

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM03451



Complete, life size leg from hip to toe, with explosion injury (blend grenade) on the lower leg and a shot right through the upper leg.

Main advantages:

- Both wounds do have a pumping mechanism and can be used for tamponade treatment.
- On the upper leg additionally training of tourniquet use is possible.
- Both wounds can bleed independently so the situation can be adjusted to your training needs.

Wound moulage "Cut"

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02820



This moulage shows a cut of about 1cm depth. The moulage adheres on human skin or a manikin, however, for longer use or for improved adherence additional skin glue is available. Supplied in storage box. Size of the moulage: approx. 12 x 8 x 1.3 cm

Wound moulage "cut with glass fragment"

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02819



This Moulage shows a cut caused by a glass fragment while the fragment is still inside the wound. The glass fragment looks realistic, but is made of soft material to avoid any risk of injury. The moulage adheres to skin or a manikin without glue, however, for longer training or improved adhesion a skin adhesive (optional) can be used. Supplied with storage box. Size of the moulage: approx. 11 x 7 x 1 cm

Wound moulage laceration, large with bleeding function

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02831



This moulage shows large lacerations. With line for bleeding function. The moulage adheres to skin or a manikin. For longer training or improved adherence an optional glue is available.

Supplied with storage box. **Size of the moulage:** approx. 28.5 x 15 x 2.3 cm

Wound moulage cut, large with bleeding function

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02830



This wound moulage shows a large, deep cut. With line for bleeding function. The moulage adheres to skin or a manikin. For longer training or improved adherence an optional glue is available.

Supplied with storage box.

Size of the moulage: approx. 28 x 14.5 x 1.8 cm

Wound moulage cut broken bottle

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02826



This moulage shows a cut by a broken bottle. The moulage adheres to skin or a manikin. For longer training or improved adherence an optional glue is available.

Supplied with storage box. Size of the moulage: approx. 10.5 x 7.2 x 1 cm

Wound moulage bullet entrance hole

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02827



This wound moulage shows an entrance hole of a bullet. The moulage adheres to skin or a manikin. For longer training or improved adherence optional glue is available.

Supplied with storage box. Size of the moulage: approx. 11 x 7.5 x 0.8 cm

Wound moulage bullet exit hole

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02828



This wound moulage shows an exit hole of a bullet. The moulage adheres to skin or a manikin. For longer training or improved adherence optional glue is available.

Supplied with storage box. **Size of the moulage:** approx. 10.8 x 7 x 1.2 cm

Wound moulage bullet hole, leg with bleeding function

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02835



This moulage shows a bullet entrance and exit wound. The moulage can be wrapped around the leg to match the two holes. With line for bleeding function. The moulage adheres to skin or a manikin. For longer training or improved adherence optional glue is available.

Supplied in storage box. Size of the moulage: approx. 29 x 14.5 x 1.5 cm

Wound moulage bullet wound, arm with bleeding function

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02834



This moulage shows a bullet entrance and exit wound. The moulage can be wrapped around the arm to match the two holes. With line for bleeding function. The moulage adheres to skin or a manikin. For longer training or improved adherence optional glue is available.

Supplied in storage box. Size of the moulage: approx. 28.5 x 15 x 1.5 cm

Wound moulage "laceration"

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02821



This moulage shows a laceration with an uneven wound edge. The moulage is self adherent on skin and on manikins. However, for longer training or improved adherence an optional glue is available.

Supplied with storage box. Size of the moulage: approx. 11.5 x 7.8 x 1.1 cm

Wound moulage "human bite"

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02823




This wound moulage shows a bite by a human. The moulage adheres to skin or a manikin, however, for longer training or improved adherence an optional glue is available.

Supplied with storage box. Size of the moulage: approx. 11 x 8 x 1 cm

Wound moulage protruding small intestines

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02841



This moulage shows protruding small intestines after injury of the abdominal wall . The moulage adheres to skin or a manikin. For longer training or improved adherence optional glue is available.

Supplied in storage box. Size of the moulage: approx. 28.5 x 14.7 x 3.5 cm

Wound moulage protrusion of large intestines

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02839



This moulage shows protruding large intestines after injury of the abdominal wall . The moulage adheres to skin or a manikin. For longer training or improved adherence optional glue is available.

Supplied in storage box. Size of the moulage: approx. 28.5 x 14.7 x 4.3 cm

Wound model - scratch

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02672



This model shows a scratch / abrasion. It is possible to add impurities (sand, small stones) and blood.

The model adheres to simulated patients or mannequins, but can also be attached to the skin with special glue to ensure longer exercise. Delivery incl. Storage box. Wound size: approx. 11 x 7.5 x 0.7 cm

Wound moulage "open fracture"

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02822



This wound moulage shows an open fracture of a smaller limb bone, e.g. lower arm. The moulage adheres to skin or a manikin, however, for longer training or improved adherence an optional skin glue is available.

Supplied with storage box. Size of the moulage: approx. 12.5 x 7.6 x 1.7 cm

Wound moulage open fracture, leg with bleeding function

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02833



This moulage shows an open fracture of a bigger limb bone (e.g. femur). With line for bleeding function. The moulage adheres to skin or a manikin. For longer training or improved adherence optional glue is available.

Supplied in storage box. Size of the moulage: approx. 28.5 x 14.7 x 3.2 cm

Wound moulage open fracture, arm bone with bleeding function

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02832



This moulage shows an open fracture of a smaller limb bone. With line for bleeding function. The moulage adheres to skin or a manikin. For longer training or improved adherence an optional glue is available.

Supplied with storage box. Size of the moulage: approx. 28.5 x 15 x 3 cm