

Low frequency acoustic pressure wave generator - SIRIO GO

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: AD03853



FEATURES • An acoustic wave with a frequency of 38 kHz • Power up to 3 W/cm² • Simultaneous connection of 2 impulse heads (flat and concave) • Possibility of cooperation with shockwave impulse head and continuous flat head (option) • Therapy of acute and chronic inflammations, calcifications and post-traumatic therapy of the musculoskeletal system • Fast and durable regeneration of soft tissues • Safe and painless therapy • Results noticeable after the first treatment

TECHNICAL PARAMETERS

- SIRIO GO
- Portable camera type
- 5.5" touchscreen screen
- Frequency 38kHz +/- 2kHz
- Maximum power 3W/cm²
- Maximum absorbed power 60 W
- Treatment time 1-30 minutes
- Device weight 2.4 kg +/- 0.2 kg
- Power supply 100/240V, 50/60Hz
- Dimensions 300x205x80mm
- Basic accessories flat impulse head Ø 50 mm (19.6 cm²),
concave impulse head Ø 50 mm (19.6 cm²)
- Optional accessories continuous flat head, shockwave head

Low frequency acoustic pressure wave generator - SIRIO

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: AD00891



SIRIO is a unique device for physiotherapy generating acoustic waves at a frequency of 38 kHz..

This revolutionary therapy works on the principle of connections of cavitation and thermal effects. This combination of nociceptive metabolites are removed from the body, which accelerates the process of regeneration and had a positive therapeutic effects.

Usage:

Sirio ideally suited for:

- Muscular Traumatology
- Muscular elongations and contractures
- Calcific Tendinitis - Duplay's disease
- Tendinopathies
- Ankle sprains
- Patellar peritendinitis
- Myositis ossificans
- Gonarthrosis
- Impingement syndrome
- Bursitis
- Low back pain
- Lateral epicondylitis (tennis elbow)
- Medial epicondylitis (golfer's elbow)
- Calcaneal spur
- Plantar talalgia
- Muscle drainage and cool-down
- Edema and hematoma

POSITIVE INFLUANCE ON tissue. SIRIO synergistically combines:

- cavitation
- THERMAL EFFECT
- MECHANICAL VIBRATIONS

THERAPY

SIRIO therapy is safe and painless for the patient. Advanced technology used in SIRIO makes therapeutic effects



complement enhancing positive action. A few minutes' treatment affects to the musculoskeletal system by increasing local metabolism and helping to remove its products.

This effect gives a quick and lasting regeneration of soft tissues. SIRIO therapy is recommended for both chronic and acute inflammatory, spurs, injuries and other diseases associated with locomotor system.

Standard equipment

- Flat silver hand probe (pulsed emission)
- Concave silver hand probe (pulsed emission)

Additional equipment

- Flat silver hand probe (continuous emission)
- Shockwave hand probe (pulsed emission)