

## Calf Disbudding Kit

**Price inquiry:** +48 605999769, [kontakt@openmedis.pl](mailto:kontakt@openmedis.pl)

Product code: AM02574



The trainer allows you to learn how to perform disbudding, remove horn buds and learn local anesthesia during the above-mentioned procedures. It is possible to learn the correct position of the nerve and inject the drug at the right angle and to the right depth. When a nerve is touched with a needle, a buzzer sounds and the LED light changes color. Once you have mastered the blockades, you can practice decorating the calf. The head of the calf has a realistic look with a hard skull, soft neck and soft reinforced ears.