

Model for learning intramuscular injections and collecting samples from the ear veins of a pig

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: MZ02663



When the correct puncture point is reached, the buzzer sounds and the LED turns green. There is also an option to turn off the buzzer during training, when the student feels they are in the right injection site, they can turn it on to see if they are right.

This anatomically correct pig ear model has a tactile silicone lateral ear vein that runs through the ear and connects to the blood pump. The ear can be folded back to expose the vein for insertion of the needle. The model can be used as one of the modules for the overall model (MZ02661)