

Wound Moulage Colostoma under skin level, incl. stand

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02966



This moulage shows a colostomy under skin level, the most common type of stoma creation. The terminal colostoma is usually found in the middle of the body or on the left side of the body and shows the colon sewn in above skin level.

Additional information:

- Due to the use of a special plastic compound, all common stoma bags adhere to the moulage.
- The model is suitable for sizing using the stoma template, cleaning the skin with water and attaching the base plate and stoma bag.
- The moulage is in a holder and can also be removed.