

Pupil simulator

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM03431



The simulator allows the user to examine the pupil reflex of the patient. Various pathologies can be simulated, for example, left or right RAPD (Relative Afferent Pupillary Defect). The simulator is easy to use, with a simple switch the reflex of the left and/or right pupil can be turned on or off.

Pupil settings:

- Both pupils constrict
- Right pupil constrict, left pupil does not react
- Right pupil does not react, left pupil constrict
- Both pupils do not react.

Technical specifications:

• 30 x 30 x 31 cm