

## CarePump Move8, 8-chamber lymphatic drainage (compression massage) device

**Price inquiry:** +48 605999769, [kontakt@openmedis.pl](mailto:kontakt@openmedis.pl)

Product code: SR03504



A modern, mobile device for lymphatic drainage

Battery powered, 5 operating modes, pressure up to 250 mmHg! The CarePump Move8 8-chamber lymphatic drainage device was designed for:

- athletes who want to achieve better results and require quick recovery before and after training;
- private users who appreciate the convenience of performing lymphatic drainage at home;
- mobile physiotherapists and those who mainly deal with sports physiotherapy;
- cosmetologists who use therapies in their offices that improve the quality of skin tissue and subcutaneous tissues and support the treatment of cellulite.

Move8 is a device that can be successfully equipped in a rehabilitation or beauty salon , but due to the battery power supply, its purpose is much wider.

Our CarePump series also includes 4-chamber and 6-chamber devices. The more chambers the apparatus supports, the more precise the lymphatic drainage treatment is, but at the same time the device is more expensive.

### Cuffs for CarePump Move8

The CarePump Move8 comes with 2 lower limb cuffs. However, you can buy other cuffs that will work with the camera, intended for upper limbs, stomach and hips, as well as a cuff in the form of shorts, which we recommend primarily for body shaping, cellulite reduction and skin firming in the area of ??thighs, buttocks and hips .

- 8-chamber wrist cuff
- 8-chamber hip and belly cuff
- Cuff - short pants 8-chamber
- Cuff - long pants 8-chamber
- 8-chamber double cuff for arms, shoulders and chest
- Cuff for the right hand, shoulder and chest 8-chamber
- Cuff for the left hand, shoulder and chest 8-chamber

All CarePump devices support "fishbone" cuffs - the "overlapping" used in them (a system of overlapping chambers) prevents the backflow of fluids during the procedure, which makes it more effective.



#### Move8 camera - operating modes

The Move8 lymphatic drainage device supports 8-chamber cuffs and has as many as 5 operating modes .

The variety of operating modes prevents the tissue from getting used to the same stimulus, thanks to which treatments performed with the device do not lose their effectiveness with prolonged and intensive use.

The device is equipped with the option of switching off individual chambers , which is especially important in the case of some diseases/injuries.

#### Kit Contents:

- Lymphatic drainage device - 8-chamber CarePump Move8
- Cuffs for the lower limbs in a universal size (in black) - 2 pcs.
- Double air duct - 1 set
- Air line distributor
- Power cable
- Convenient transport bag.
- Features of the device
- The highest quality
- The highest quality 8 chambers in the cuffs
- 8 chambers in the cuffs 5 operating modes
- 5 operating modes Compact design
- Compact design Pressure regulation in the range of 30-250 mmHg
- Pressure regulation in the range of 30-250 mmHg Possibility to set the treatment time from 10-60 minutes
- Possibility to set the treatment time from 10-60 minutes Battery operated
- Battery powered Possibility to turn off individual chambers
- Possibility to turn off individual chambers Low volume
- Low volume

The treatment primarily facilitates the removal of excess fluid from the muscles, including the removal of metabolic metabolites, as well as inflammatory mediators (due to loads). Such action leads to relaxation of the tissues and reduction of possible ailments.

Lymphatic drainage supports the treatment of lymphedema - it serves to "push" the accumulated fluid and eliminate swelling by improving lymph circulation. The purpose of the massage is both to prevent the formation of stagnation and to help remove the resulting ones.

#### Indications for presotherapy treatments (lymphatic drainage):

- pre- and post-workout regeneration,
- prevention of circulatory and lymphatic system failure,
- diseases causing swelling,
- feeling of "heavy legs",
- cellulite,
- exudates,
- prevention and prophylaxis in the treatment of obesity and overweight,
- past infectious diseases (in order to strengthen the body and remove toxins),
- accelerating the absorption of hematomas or bruises by stimulating microcirculation.

#### Lymphatic drainage in the regeneration of athletes

Lymphatic drainage devices support the flow of blood and lymph, thanks to which they significantly accelerate regeneration. For this reason, they are willingly used by athletes before training or competing in order to achieve better results, as well as after exercise in order to regenerate (renew) faster and maximize training effects.

Lymphatic drainage is an indispensable element of regeneration at training camps and during the peak of the sports season, when the player has to keep fit for a longer period of time.

The Move4, Move6 and Move8 compression massage devices are ideal for athletes due to the fact that they have as many as 5 different operating modes and a large range of pressure regulation generated by the device. Battery



power allows you to perform the procedure anywhere, also outdoors or when returning from competitions.

Recommendations for working on the battery:

The battery included in the equipment of the device is a consumable element - the number of battery charges is not less than 500. We understand the charging cycle by connecting the charger to the device each time. To prolong the life of the battery, it is recommended to discharge it until the red indicator light comes on. Then complete the full charging cycle again.

Dimensions of the cuffs for the lower limbs

- Thigh circumference (at the widest point) - 70 cm
- Thigh circumference (center) - 59 cm
- Knee circumference - 52 cm
- Calf circumference - 48 cm
- Ankle circumference - 44 cm
- Thigh - heel - 87 cm
- Thigh - metatarsal - 102.5 cm
- Toes - heel - 22.5 cm