

Cervical Dilatation & Effacement Module (CD&E) - PROMPT Flex

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM03664



Part of the PROMPT Flex range, the Cervical Dilatation & Effacement Module allows for training in the assessment and management of both the latent, and active, early stages of labor. This model can be used to create a wide variety of birthing scenarios to identify:

- Position/tilt of the cervix
- Consistency of the cervix
- Cervical effacement
- Cervical dilation – early/latent phase, active phase and transition phase
- Fetal station
- Normal and abnormal foetal presentation
- Moulding of the skull bones
- Caput formation
- Presence of membrane and amniotic fluid

This cervical dilation module can be purchased as an add-on to the Standard or Advanced PROMPT Flex simulators, and allows for the assessment of labor progression, likelihood of vaginal delivery, and suitability for induction. **SKILLS GAINED**

- Vaginal examination to assess:

Position of the cervix
Foetal presenting part
Moulding and caput formation
Presence/absence of membrane

- Calculation of Bishop's Score (also known as Pelvic Score):
- 5 measurements of the pelvic examination:

Centimetres dilated
Effacement
Foetal station
Consistency
Position

- Other score modifiers:

Pre-eclampsia
Each previous vaginal delivery
Postdate pregnancy
Nulliparity
Overweight
Preterm rupture of membranes

- Trainees are able to practise the skills required to assess both normal and complicated births



- Palpable ischial spines allowing for accurate assessment of station
- Model accommodates a variety of training methods

REALISM

- Rupturable membrane, with fluid, on flexed or deflexed head
- Dynamic Rig allows adjustment of dilation enabling scenario based, and communication skills, training
- Realistic representation of cervixes, including anterior lip

VERSATILITY

- Offers a large variety of presentations, including: stations, relative positions, and degrees of rotation
- Easy removal and replacement of presenting parts and cervixes when using the Static Rig
- Presenting parts and cervix can be set up externally to the model for demonstration

CLEANING

- Remove excess lubricant from the component surfaces
- Clean components with a soft damp cloth, warm water and mild detergent
- Allow to dry before storing trainer

SAFETY

- This product is latex free
- Using alternative lubricants to that recommended, may weaken or damage components of the cervical dilation model
- Wear gloves while handling the simulator and during procedures, sharp surfaces (fingernails and jewellery) can lead to silicone damage
- To prolong the life of caputs, keep covered in the supplied plastic tray when not in use

ANATOMY

- Cervices:

0cm dilation, unripe/firm, no effacement

1cm dilation, unripe, 40% effacement

2cm dilation, medium ripeness/firmness, 50% effacement

3cm dilation, ripe/soft, 65% effacement

4cm dilation, ripe/soft, 70% effacement

4cm dilation, ripe/soft, 100% effacement

Anterior lip, approximately 10cm dilation

Dynamic cervix for use with dynamic rig, allows for cervical dilation from 4cm to 10cm, ripe/soft, 100% effacement

- Presenting parts:

Flexed head

Deflexed head

Face presentation

Brow presentation

Breech presentation

- Moulding and caput formation:

Flexed or deflexed head

Parietal bones, one fixed and two moveable to simulate variation of moulding

Caput

- Perineum for CSM and CD&E Modules - PROMPT Flex
- Pelvic Ring for CD&E Module - PROMPT Flex
- Latent Labor Cervixes for CD&E Module - PROMPT Flex
- Dynamic Cervix for CD&E Module - PROMPT Flex
- Presenting Parts for CD&E Module - PROMPT Flex



- Difficult Labor Heads for CD&E Module - PROMPT Flex
- Amniotic Membrane Set for CD&E Module - PROMPT Flex
- Caputs for CD&E Module - PROMPT Flex

WORKS WITH

- Birthing Simulator PROMPT Flex - Standard (Light Skin Tone)
- Birthing Simulator PROMPT Flex - Advanced (Light Skin Tone)