


Forelimb Canine Bandaging & IV Kit

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: MZ03805



An IV and Bandaging training model all in one.

Veterinary students can practise finding the cephalic vein in a controlled environment for intravenous therapy or blood sampling with clear assessment and feedback from the tutor. Blood circulates from the Holsim Blood Pump through the vein tubing, where the student can then attempt to identify the cephalic vein and withdraw blood or administer the required fluids/drugs. The leg can also be used for teaching bandaging methods, especially as the toes of the dog leg can be spread to wrap the bandage. Forelimb Canine Bandaging & IV Kit Includes:

- 1 x IV & Bandaging Forelimb Model with IV Insert
- 1 x Blood Pump
- 1 x Jar of Artificial Blood

About

Anatomically Correct

- Bone and vein positioning have been replicated from real canine anatomy.

Blood Pump

- Blood is circulated through the leg by a small electronic pump to simulate the heart pumping.

Needle Resistance

- When inserting the needle, a realistic resistance is felt under the skin.

Bulging Vein

- When pressure is applied, the vein bulges for students to see where to insert the needle.

Spreadable Toes

- Toes can be spread for bandaging purposes between digits.