

IV Insert - Hindlimb Bandaging Model

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: MZ03808



Designed to be used with the Holsim Hind Limb Bandaging & IV Model.

The dog skin insert is first connected to the tubes in the dog leg. The blood will circulate through the vein tubing, where the student can then attempt to identify the vein and withdraw blood or administer the required fluids/drugs. The insert has a realistic resistance under the skin when the needle is inserted. The vein bulges slightly when pressure is applied. The insert also has coloured skin, fat and muscle layers.