

K9 Otoscopy Simulator

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: MZ04029



Our otoscopy simulator can be used to train how to correctly carry out an ear examination in dogs. The external auditory canal up to the eardrum is realistically reproduced. The ear canals were obtained from CT data. A healthy as well as a narrowed ear canal is modeled. To check success, the eardrums in both ear canals were marked in color. There is an opening in the healthy ear canal so that artificial cerumen can be introduced into the model. The dog head is made entirely of soft silicone and is detachably mounted on a tripod. The tripod offers the possibility of fixing the dog's head in a frontal and side examination position. The table clamp allows the tripod to be firmly attached to a table.