

Injection Belly

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM01465




This soft stomach replica will be an invaluable tool for teaching patients how to properly insert and rotate the infusion set for insulin pumps, avoiding the 5 cm area around the navel, as well as giving self-injections of many types including pen-style injections. The trainer has lifelike skin with what feels like real stomach tissue underneath. There is a soft plastic backing to keep needles from going through and holes to allow the material to breathe and dry, should you choose to inject liquids. Distilled water may be injected, but the more liquid that is injected into the simulator, the longer the time for evaporation. Actual insulin or other medications should not be used, as they may deteriorate the material. The injection belly is considered disposable; however, with proper care and treatment it will serve its purpose for a long period of time. Using the smallest needles possible when simulating injections and rotating the insertion points will help prolong the life of the trainer.

Size: 25.4 x 15.3 x 3.2 cm, Weight: 0.4 kg

Diabetic injection pad, advanced version

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM01854



Unlike other injection pads this advanced pad gives the opportunity to practice with a normal and a hardened abdominal wall. It allows the student to form a realistic fold for injection. Due to special material it gives a lifelike impression of touch and feel, almost like in a real patient. The model is perforation proof and can be attached to the students' body using a flexible band.

Advanced Injection Belly

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM01853



This soft stomach replica is an invaluable tool for showing how to properly insert and rotate the infusion set for insulin pumps, avoiding the 5 cm area around the navel, as well as giving many types of self-injections (including pen-style injections). Similar to the Injection Belly, but made of a material that creates a more realistic experience. The trainer has lifelike skin with what feels like real stomach tissue, along with the ability to pinch the skin. Backing keeps needles from going through. Adjustable straps allow the stomach to be attached to a manikin or human patient for simulation practice. The soft material can be punctured repeatedly without showing holes. Liquid injection is not recommended.

Size: 25.4 x 15.3 x 5.1 cm, Weight: 0.4 kg