

Median section of a dog head

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: M01293



This scaled down model of a median sectioned dog head depicts graphically the anatomy of mouth, tongue as well as throat with trachea and esophagus. Also the olfactory organ, the skull with sinuses and a brain section with cerebellum, brainstem and afterbrain are visible. Supplied on stand.

Size without stand: 17 x 15 x 6 cm, Weight: 0.6 kg

Canine Jaw healthy/diseased

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: M01358



This canine jaw of a mid-size dog shows healthy teeth on the right side and diseased/broken teeth on the left side. The nine pathologies shown in this model are: missing premolar, gingival recession, retained deciduous tooth, worn incisors, gingivitis, plaque, tarter accumulation, periodontal disease and fractured canine. The jaw can be opened, closed and separated for detailed study.

Size: 12 x 7.5 x 5.5 cm

Canine Ear - Normal/Infected

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: M01360



This 2-sided, average size, canine ear depicts: a normal side with cochlea, auditory ossicles, auditory tube, tympanic bulla, middle ear cavity, tympanic membrane, horizontal canal, vertical canal, auricular cartilage, pinna and temporalis muscle; abnormal side illustrates inflamed inner ear structures, inflammatory exudate in tympanic bulla, ear canal with partial occlusion from cellular hyperplasia, inflammatory exudate and an inflamed (reddened) outer ear.

Size: 12 x 7 x 5 cm

Canine Parasites

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: M01402



It is fairly common for a dog to become infected with an external parasite at some point in its lifetime. Parasites can affect your dog in a variety of ways, ranging from simple irritation to causing life-threatening conditions if left untreated. Some parasites can even infect and transmit diseases to people. This model here shows in great detail the three most common external parasites a mite, flea, and tick.

Size: appr. 27 x 16 x 15 cm

K9 Otoscopy Simulator

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: MZ04029



Our otoscopy simulator can be used to train how to correctly carry out an ear examination in dogs. The external auditory canal up to the eardrum is realistically reproduced. The ear canals were obtained from CT data. A healthy as well as a narrowed ear canal is modeled. To check success, the eardrums in both ear canals were marked in color. There is an opening in the healthy ear canal so that artificial cerumen can be introduced into the model. The dog head is made entirely of soft silicone and is detachably mounted on a tripod. The tripod offers the possibility of fixing the dog's head in a frontal and side examination position. The table clamp allows the tripod to be firmly attached to a table.