


## Retracted Ostomy Care Training Model

**Price inquiry:** +48 605999769, [kontakt@openmedis.pl](mailto:kontakt@openmedis.pl)

Product code: SM01988



Designed to help introduce the essentials of ostomy care to patients and students, this model is ideal for demonstration and practice of ostomy care procedures on a stoma that has sunk below skin level, including cleaning and changing the pouching system. The anatomy of the ostomy is carefully reproduced to provide lifelike functions and appearance. A soft, pliable material is used for the stoma in order to achieve the most realistic tactile sensation. Dilation of the stomas can be demonstrated and practiced, along with application of postoperative and permanent ostomy bags on stomas that can often be difficult to obtain a good seal around. The ostomy can actually be irrigated. Syringe

is included to pump simulated stool through the unit to provide drainage and excretion at the ostomy. The consistency of the stool can be varied by using water to thin the material.

Includes model, display easel, syringe, simulated stool, lubricant, and instruction manual.

Size: 18 x 8 cm, Weight 1.6 kg