


Nasogastric Equine Intubation Model

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: MZ02557



A realistic trainer designed for learning and training nasogastric intubation in horses. Anatomically correct structure reflects realistic treatment sensations. It allows you to perform basic surgery training routinely performed in equine colic cases to decompress the stomach and provide therapy. This model can be secured with the head facing down first to avoid inserting the tube into the trachea. The nasogastric tube can then be inserted into the nose in the ventral and medial directions. As the tube reaches the throat, there is soft resistance as the tube contacts the esophageal opening. The horse's head can be bent up and down to open the airway or esophagus. The model includes both the esophagus and the trachea to ensure proper tube placement. You can touch the nasogastric tube above the left jugular vein and feel it travel along the neck under the skin to the stomach. The model also includes dorsal ("false nostrils") and abdominal ("true nostrils").

Attention!

The photos are illustrative and the product may finally differ in appearance, e.g. colors.

The skull is produced in white or gray depending on the availability of materials.