

Human Vascular Access Sleeve Light Skin

Price inquiry: +48 605999769, kontakt@openmedis.pl

Product code: SM02958



Train on venipuncture and vascular access skills in life-like scenarios with no harm to participants with the Human Vascular Access Sleeve. Designed to fit most people and manikins, the sleeve is secured by the hook and loop tabs above and below SurgiReal's simulated RealFlow tissue pad.

Keeping safety as our top priority, all Vascular Access Sleeves are designed with puncture resistant material to make sure the needle won't make contact with your participant.

Training applications include, but are not limited to:

- Needle Alignment and Vascular Access,
- Palpating the Vessel,
- IV Placement.

BUILD YOUR SKILLS IN THREE STEPS:

Sleeve attached to manikin,
Sleeve attached to stationary person,
Sleeve attached to a person mimicking real-life patient movement.

- 1 Sleeve,
- 1 Attachable Simulated Tissue Pad,
- 1 Package of Artificial Blood Powder,
- 1 IV Bag (1000 ml).