

## Small leg support with skin crack

**Price inquiry:** +48 605999769, [kontakt@openmedis.pl](mailto:kontakt@openmedis.pl)

Product code: MZ03786



Train stitching lacerations and stitching on a living object. The leg mounts are suitable for most small animals and training dummies and are secured by the Velcro fasteners above and below SurgiReal's simulated RealLayer Tissue Pad. Possible applications include practicing superficial seams, ensuring knot safety and sinking seams. Since safety is our top priority, all skin crack holders are made of a stab-proof material, so that the needle cannot come into contact with your participant. Develop your skills in three stages: holder attached to a manikin, bracket attached to a stationary actor, bracket attached to an actor imitating the movement of a patient